


CRYOGENIC RECEIVERS



Part No: G - 050 - E - LN6N

Operating Note

This unit is a high performance cryogenically operated germanium photodiode/amplifier designed for low frequency DC or chopped measurements. The output voltage is proportional to radiation incident on the active area as follows:

 $V_{out} = P_{sig} \times R_{\lambda} \times R_{f}$

where P_{sig} is incident power in watts, R_{\(\lambda\)} is the photodiode responsivity in A/W at the wavelength of interest, and R_f is the amplifier transimpedence gain. This is DC coupled with high gain and extensive care should be taken in shielding the unit from any ambient light during operation. Exposure to room lights may cause amplifier saturation and can the lead to failure of unit.

SPECIFICATIONS	
Active Area	5 mm diameter
Spectral Range	800 – 1800 nm @ 298K; 800-1500 nm @ 77K
Shunt Resistance	> 10 kΩ @ 298K; > 1000 MΩ @ 77K
Shunt Capacitance	30000 pF typical
NEP	$< 1.5 \times 10^{-12} \text{ W/Hz}^{1/2} @ 298\text{K}; < 1.5 \times 10^{-15} \text{ W/Hz}^{1/2} @ 77\text{K}$
Responsivity @ 1.3 μm	0.8 A/W minimum, 0.9 A/W typical
Dewar Hold Time	12 hours minimum with liquid N_2
Field of View	60° nominal
Amplifier	Dual-Gain Transimpedance
Gain	1 x 10 ¹⁰ V/A / 1 x 10 ⁹ V/A
Bandwidth	DC - 15 Hz / DC - 150 Hz
Connections	BNC signal coaxial cable with 3 lead shielded power cable. Red = +V, Black = -V, White/Shield = ground Note: A DB9 connector is provided on units purchased with optional PS-1 Low Noise Power Supply